



Revision 5:
GeNUGate 600/800 (ALG), GeNUCenter: 600/800,
GeNUScreen: 600/700/800, GeNUBox: 600/800

em							
0	1						

em		0					
2	3	1					

em		0	2				
4	5	1	3				

em		0	2	4			
6	7	1	3	5			

em		0	2	4	6		
8	9	1	3	5	7		

em		0	2	4	6	10	
8	9	1	3	5	7	11	

em		0	2	4	6	10	12
8	9	1	3	5	7	11	13